



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fixed Bed Adsorption Plant

Product Code :
LBNY-0005-11600026



Description :

Adsorption column with hot water jacket for temperature control
Hot water generator with pumps and electronic thermostat
3x electro-pneumatic converters
2 Thermal mass flow meters with 4
3x pneumatic manual control valves
20 mA output signal
6x digital thermocouples with pt100 probes protected by stainless steel shield sleeves
Infrared transmitter
Pressure transmitter
Electrical control panel in main cabinet with earth leakage circuit breaker
Pipework manufactured from stainless steel
Microprocessor digital PID controller with serial card
Emergency stop button
Supervision and control software for Windows.

Technical Specification :

Fixed Bed Adsorption Plant

Laboratory Instrument India

