



Email : sales@laboratoryinstrumentindia.com

Product Name :
Manual Tray Dryer Process Unit

Product Code :
LBNY-0005-11600027



Description :

Horizontal air duct with 450 mm x 450 mm window.
Axial fan with digital speed control 0 to 4 meters per second.
2x digital thermohygrometers with digital displays.
4x aluminium drying trays 400 mm x 300 mm.
Electrical heating element array 2,700 watts.
Distillation column height 600 mm with 7 sieve trays.
Electronic balance with range 0 to 10 Kgs and 0.1 g resolution.
Digital anemometer with digital display.
Digital chronometer.
Emergency stop button
Electrical control panel in main cabinet.

Technical Specification :

Manual Tray Dryer Process Unit



Laboratory Instrument India

