



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Manual Tray Dryer Process Unit

**Product Code :**  
LBNY-0005-11600027



**Description :**

Horizontal air duct with 450 mm x 450 mm window.  
Axial fan with digital speed control 0 to 4 meters per second.  
2x digital thermohygrometers with digital displays.  
4x aluminium drying trays 400 mm x 300 mm.  
Electrical heating element array 2,700 watts.  
Distillation column height 600 mm with 7 sieve trays.  
Electronic balance with range 0 to 10 Kgs and 0.1 g resolution.  
Digital anemometer with digital display.  
Digital chronometer.  
Emergency stop button  
Electrical control panel in main cabinet.

**Technical Specification :**

Manual Tray Dryer Process Unit



**Laboratory Instrument India**

