



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spray Dryer Process Unit

Product Code :
LBNY-0005-11600028



Description :

Temperature programmable up to 2500C
Touch screen control panel
Variable speed fan
Built-in compressor
Heating system 3,000 watts
Atomization system with spray nozzle and cleaning device
Peristaltic pump adjustable up to 1,500 ml per hour
Exhaust pipe for solvent vapour
Cyclone with bottle for collecting the dried product
Drying chamber with collection vessel for residue

Technical Specification :

Power supply: 230 Vac 50 Hz single-phase -3000 VA.



Laboratory Instrument India

