



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Spray Dryer Process Unit

**Product Code :**  
LBNY-0005-11600028



**Description :**

Temperature programmable up to 2500C  
Touch screen control panel  
Variable speed fan  
Built-in compressor  
Heating system 3,000 watts  
Atomization system with spray nozzle and cleaning device  
Peristaltic pump adjustable up to 1,500 ml per hour  
Exhaust pipe for solvent vapour  
Cyclone with bottle for collecting the dried product  
Drying chamber with collection vessel for residue

**Technical Specification :**

Power supply: 230 Vac 50 Hz single-phase -3000 VA.



**Laboratory Instrument India**

