



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hemoglobino-Meter

**Product Code :**  
LBNY-0001-510004



**Description :**

Hemoglobino-Meter

**Technical Specification :**

The portable hemoglobin meter is a micro processed photometer. In a disposable vial, containing Drabkin's reagent, 10 uL of blood sample are dropped .

Hemoglobin content is read and calculated by a microprocessor and proprietary software. Results are presented in a LCD display.

Sample type: capillary blood & venous blood

Wavelength: Dual wavelength, one for hemoglobin measurement and one for turbidity compensation

Size and Display system: Compact and lightweight portable Hemoglobinometer with Digital display

Display range: Measure hemoglobin in the range of 0-30 g/dl

Measuring time: measure hemoglobin within 1 minute

Sample volume: Minimum sample volume up to 50 µl

Result storage: Minimum of 500 results with date and time, have data transfer facility to printer or computer (optional )

Printer: Built in Printer (optional)

Operating Temperature: +15°C to +30°C)

Utility Requirement: Electrical power supply: 220VAC ±10%, also battery operated

Hemoglobino-Meter, Hemoglobino-Meter Equipments, Hemoglobino-Meter Tools, Hemoglobino-Meter Tool Kits, Hemoglobino-Meter Manufacturers, Hemoglobino-Meter Suppliers from India, China, Kenya



**Laboratory Instrument India**