



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Domestic Conditioning Trainer

**Product Code :**  
LBNY-0005-11700011



**Description :**

Domestic Conditioning Trainer allows to study the split type conditioners with inverter. The system is supplied complete with technical manual and experiments. It enables the student to acquire a good knowledge and practice with the main components of this kind of equipment.

**Technical Specification :**

Electrical supply: 220 V, 50/60 Hz, single-phase, 2 Kw

Dimensions and Weight:

Dimensions: 1600x800x1900 h mm

Weight: 100 kg.



**Laboratory Instrument India**