



Email : sales@laboratoryinstrumentindia.com

Product Name :
Domestic Conditioning Trainer

Product Code :
LBNY-0005-11700011



Description :

Domestic Conditioning Trainer allows to study the split type conditioners with inverter. The system is supplied complete with technical manual and experiments. It enables the student to acquire a good knowledge and practice with the main components of this kind of equipment.

Technical Specification :

Electrical supply: 220 V, 50/60 Hz, single-phase, 2 Kw

Dimensions and Weight:

Dimensions: 1600x800x1900 h mm

Weight: 100 kg.



Laboratory Instrument India