



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vacuum And Recharge Station

**Product Code :**  
LBNY-0005-11700012



**Description :**

The Vacuum And Recharge Station unit allows the assembling and disassembly of refrigeration and air conditioning plants and the filling and draining operations of the refrigerating fluid.

Vacuum and Recharge Station Composed of:

Portable recharge station with trolley-mounted frame, fitted with two-stage pump 120W for high vacuum, graduated charging cylinder

The two-stage pump used is fitted with a ballast valve and anti-air return valve

a. Flow rate: 3.6 m<sup>3</sup>/h

b. Cylinder capacity: 1.1 kg, capacity more than enough for filling of our educational equipment

c. Dimensions: 450 x 280 x 550 h mm

Ratchet wrench for refrigeration

Multigas electronic leak detector able to rapidly locate gas leaks.

Adjustable wrench max. Opening 30 mm

Cylinder connector.

**Technical Specification :**

Vacuum And Recharge Station

---

# Laboratory Instrument India

