



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Set Of Sectioned Components For Refrigeration And  
Conditioning Plants

**Product Code :**  
LBNY-0005-11700015



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The Set Of Sectioned Components For Refrigeration And Conditioning Plants set includes some of the most important components, used in refrigeration and conditioning plants such as compressors, condensers, expansion valves etc. Special explanatory technical plates show each sectioned component explaining its operation with notes, drawings and technical data. In sectioning each component the best care has been taken to evidence all the interesting constructive detailed features. The use of contrasting colors according to a pre-established code simplifies the informative process as well as the learning one.

**Technical Specification :**

Includes:  
Two cylinder open type refrigeration compressor 1.5 HP  
Hermetic refrigeration compressor  
Absorption liquid chiller (15 l)  
Refrigeration gas condenser with concentric pipes  
Semi-hermetic refrigeration compressor 1 HP  
Air conditioning unit 1 HP.

---

# Laboratory Instrument India

