



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigeration Mobile Lab

Product Code :
LBNY-0005-11700022



www.laboratoryinstrumentindia.com

Description :

Compact laboratory designed for training of 14 years students and over. The trolley is completely mobile and has 5 storage drawers, protecting the equipment in shock-free, molded sponge pods. You can organize the equipment according to the subjects for quick and easy setup and control of the equipment itself. It is particularly suitable to introduce the main and practical refrigeration applications in both vocational and technical schools. It is possible to compose the circuit on the aluminum surface to place on the working top. Opportunity of fault insertion and detection. It is ideal for all teaching environments and its modular design guarantees both flexibility and adaptability for all the experiments. The Refrigeration Mobile Lab is already equipped with all the necessary feed fittings according to the industrial standards.

Technical Specification :

Electrical power supply: 220 V, single-phase, 50/60 Hz.

Weight and dimensions:

Weight: 300 kg

Overall dimensions: 1500x600x930 mm.



Laboratory Instrument India

