



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Refrigeration And Conditioning

**Product Code :**  
LBNY-0005-11700025



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The Refrigeration And Conditioning has been designed for both theoretical and practical learning of refrigeration technics. The unit has been designed for operation with CFC free Freon as primary and unfreezable mixture as secondary refrigerant. On the front panel, a synoptic diagram reproduces the plant clearly showing each component and the measurement points provided in the circuit. The Refrigeration And Conditioning unit is self-sufficient and provided with all the instruments which are necessary to allow a wide experimentation.

The Refrigeration And Conditioning unit is available in two versions:  
Refrigeration Study Unit with Heat Pump.  
Refrigeration Study Unit.

**Technical Specification :**

Refrigeration And Conditioning



**Laboratory Instrument India**