



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Refrigerating Room

**Product Code :**  
LBNY-0005-11700026



#### **Description :**

The Refrigerating Room is an extremely easy-to use and functional unit able to permit the study of the characteristics of a refrigerating room to store supplies at low temperature. The unit is supplied with a manual which describes the components, the installation and utilization procedures as well as many exercises with experimental results. Front panel with mimic diagram and instrumentation allows to keep the most significant thermodynamic quantities under control, so facilitating the students understanding and the teachers task.

#### **Technical Features:**

Compressor power: 150 W  
High temperature thermostat:  $-5 \text{ }^{\circ}\text{C}$  -  $+25 \text{ }^{\circ}\text{C}$   
Low temperature thermostat:  $-7 \text{ }^{\circ}\text{C}$  -  $-35 \text{ }^{\circ}\text{C}$   
Refrigerating room capacity: 30 l

#### **Technical Specification :**

Electric supply: 220V single-phase, 50Hz.  
Weight and Dimensions;  
Dimensions: 600x600x1350 h mm  
Net weight: 95 Kg.

---

# Laboratory Instrument India

