



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Air Conditioning Mobile Lab

**Product Code :**  
LBNY-0005-11700029



**Description :**

Compact laboratory designed for training of 14 years students and over. The Air Conditioning Mobile Lab is already equipped with all the necessary feed fittings according to the industrial standards. It is ideal for all teaching environments and its modular design guarantees both flexibility and adaptability for all the experiments. It is particularly suitable to introduce the main and practical air conditioning applications in both vocational and technical schools. The trolley is completely mobile and has 5 storage drawers, protecting the equipment in shock-free, molded sponge pods. It is possible to compose the circuit on the aluminum surface to place on the working top. Opportunity of fault insertion and detection. You can organize the equipment according to the subjects for quick and easy setup and control of the equipment itself.

**Technical Specification :**

Electrical power supply: 220 V, single-phase, 50/60 Hz.

Weight and dimensions:

Overall dimensions: 1500x600x930 mm

Weight: 300 kg.



**Laboratory Instrument India**

