



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Tube Roller

Product Code :
LBNY-0001-510007



Description :

Laboratory Tube Roller

Technical Specification :

Six roller design for gentle rolling and easy to clean

Continuous or timer operation.

Minimum Number of rollers: 6

Minimum Speed: 5rpm

Minimum load: 1kg

Controls type : digital

Compatible with different sizes of tubes

Timer: Required

Operating Environment: Operating temperature 15? to 30?

Laboratory Tube Roller, Laboratory Tube Roller Equipments, Laboratory Tube Roller Tools, Laboratory Tube Roller Tool Kits, Laboratory Tube Roller Manufacturers, Laboratory Tube Roller Suppliers from India, China, Kenya

Laboratory Instrument India

