



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital Laboratory Timer

**Product Code :**  
LBNY-0001-510008



**Description :**

Digital Laboratory Timer

**Technical Specification :**

Loud Sounding Alarm - You will easily be able to hear this alarm in another room. It is loud and clear but not deafening.

The mildly hard of hearing should have no problem with the alarm

Magnetic Back, Retractable Stand and Hook for Hanging - You have 3 different choices for placement of the Timer.

Either stick it to the fridge door with its large magnet back, or use its kickstand in the back for table-top use, LCD display: hours, minutes and seconds display

Timer unit: hours, minutes, seconds

Counts: Counts Up and Down timer

Alarm: Loud, Audible Alarm

Switch: On/Off Switch to increase Battery Life

Battery type : AA or AAA

Supplied with a 3-year battery

Digital Laboratory Timer, Digital Laboratory Timer Equipments, Digital Laboratory Timer Tools, Digital Laboratory Timer Tool Kits, Digital Laboratory Timer Manufacturers, Digital Laboratory Timer Suppliers from India, China, Kenya



**Laboratory Instrument India**