



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Flow Round A Bend Investigation

**Product Code :**  
LBNY-0005-11700055



**Description :**

A 90 degree bend of constant cross section has air blown through it from the Airflow System. When connected to a suitable multi-tube manometer the static pressure profile along the inner and outer radius of the bend may be measured at a range of air velocities. Along the inner and outer radius at strategic points, static pressure tapping points are located.

**Technical Specification :**

Flow Round A Bend Investigation



**Laboratory Instrument India**