



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Airflow System Base Unit

**Product Code :**  
LBNY-0005-11700056



**Description :**

A small footprint, high volume high pressure centrifugal fan with adjustable flow control, inlet and outlet couplings.

A comprehensive range of optional accessories are available that may be connected to the fan inlet, outlet or both as required.

11 Optional Heat Exchangers available:-

Flow Visualization Investigation

Multitube Manometer

Bernoulli's Equation

Flow Around a Bend Investigation

Jet Attachment Investigation

Boundary Layer Investigation

Round Turbulent Jet Investigation

Drag Force Investigation

Principles of Airflow, Friction Losses in Bends and Pipe Elements

Principles of Airflow, Fan Test and Flow Measurement

Principles of Airflow, Pressure and Velocity Distribution (Pitot-Traverse).

**Technical Specification :**

Airflow System Base Unit



**Laboratory Instrument India**