




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Internal Combustion Engine Test Stand

**Product Code :**  
LBNY-0005-11700057



**Description :**

A separate instrumentation console allows measurement and display of the engine torque, speed, air consumption, fuel consumption, inlet air temperature and exhaust gas temperature.

A floor mounted modular internal combustion engine test stand with regenerative dynamometer and a range of optional internal combustion engines up to 4kW shaft power.

The unit is supplied with a detailed experimental operating and maintenance manual giving example experimental results and sample calculations.

An optional computerized data acquisition upgrade is available that allows all relevant system parameters to be automatically recorded on a PC.

Accessories and spares for two years normal operation together with a full two year warranty are also included as standard.

**Technical Specification :**

Internal Combustion Engine Test Stand



**Laboratory Instrument India**

