



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bench Top Hot Plate for Laboratory

Product Code :
LBNY-0001-510009



Description :

Bench Top Hot Plate for Laboratory

Technical Specification :

Aluminium Moulded Construction
Highly scratch and chemical resistant Top
Even heat distribution on surface of Top
Large top plate to assist beaker upto 5000 ml
High Magnetic strength
Single multi-turn intelligent knob operation
User settable High temperature cutoff and hard wired cutoff

Technical specification:

Feature Non-slip and chemical resistant top plate and easy to clean for heating lab reagents and other purposes
Body Material Surface Type: Ceramic :or aluminum chemical and scratch resistance ceramic plate,
Temperature control: Microprocessor PID temperature control with Hot top warning indicator
temperature uniform across the plate,
Max Surface Temp °C: 540
Display: Digital temperature display with a LCD/LED indicator

Temperature Range: 50°C to 350°C
Capacity : up to 10 liters,
Dimention: L * W * H = 41.1 * 28.7 * 10.2 cm

Power supply: 220V +10% 50/60HZ

Bench Top Hot Plate for Laboratory, Bench Top Hot Plate for Laboratory Equipments, Bench Top Hot Plate for Laboratory Tools, Bench Top Hot Plate for Laboratory Tool Kits, Bench Top Hot Plate for Laboratory Manufacturers, Bench Top Hot Plate for Laboratory Suppliers from India, China, Kenya



Laboratory Instrument India