



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multi Tube Mano Meter

Product Code :
LBNY-0005-11700062



Description :

The device allows up to 16 pressures to be monitored simultaneously either relative to atmospheric pressure or another pressure, via the common reservoir. In order to increase the device sensitivity the device may be inclined at a known angle.

A multi-tube manometer with a common reservoir, that may be used to give a graphic display of pressure distribution on multi-point pressure tapings.

Technical Specification :

Multi Tube Mano Meter



Laboratory Instrument India