



Email : sales@laboratoryinstrumentindia.com

Product Name :
Environmental Chamber

Product Code :
LBNY-0005-11700065



www.laboratoryinstrumentindia.com

Description :

This precisely designed laboratory equipment is also suitable for executing cold resistance test, thermal cycle test and RH test. Test section of this chamber is designed with wind board that regulates flow direction of air for better air circulation purpose. Developed from 304 grade stainless steel, it contains platinum made temperature sensor, easy to dismantle shelves, multiple air circulating fans, touch screen based controller and LAN based network connectivity. We are an eminent manufacturer and exporter of Environmental Chamber.

Features:

Application of mechanical compression based refrigeration system
Coated outer structure
Protected against over temperature
Mechanical water supply arrangement.

Technical Specification :

Environmental Chamber



Laboratory Instrument India

