




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Zoom Flumes And Gauging Weirs Apparatus

**Product Code :**  
LBNY-0005-1180005



**Description :**

Parshall flumes constructed in GRP, supplied with measuring scales, spirit levels and calibration curves.  
Broad crested weir constructed of durable glass reinforced plastic.  
Triangular shaped broad crested Crump weir constructed of durable glass reinforced plastic.  
Stilling well with depth gauge comprising hook and point gauge.  
Thin plate weirs constructed of stainless steel and mounted on a painted weir plate carrier with a simple head scale attached.  
WSC flumes constructed in GRP and fitted with depth measuring scale supplied with calibration curve.  
Lightweight portable channel section constructed of corrosion resistant material.

**Technical Specification :**

Zoom Flumes And Gauging Weirs Apparatus



**Laboratory Instrument India**