



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drain Permea Meter

Product Code :
LBNY-0005-1180006



Description :

The Drain Permeameter is suited for use both as a teaching and demonstration tool and for laboratory testing and research. A removable base section is provided into which the filter material to be investigated is placed. The equipment consists of a 100mm diameter transparent column which may be filled with any soil type. The apparatus is supported on a sturdy metal frame and may be conveniently placed on a laboratory bench close to a water supply and drain. Water is added to the column under a constant head and permeating water and soil is collected at the base. The discharge from the drain can be monitored to assess the efficiency of the filter material.

Technical Specification :

The apparatus, which may be bench-mounted, is supplied complete with user instructions and manuals. Acrylic column, 100 mm diameter supported on a metal stand;
Removable test section is provided at the base of the column to house the filter medium to be tested.



Laboratory Instrument India