



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Water Bath

Product Code :
LBNY-0001-510011



Description :

Laboratory Water Bath

Technical Specification :

Stainless Steel, insulated double walled,
Inner wall of stainless steel,
Aluminium /SS cover,
Brass drain cock,

Temperatures setting: Digital microprocessor display to set temperature point preventing thermal
Thermostatic temp. control from ambient to 85 -100° C (±0.50) complete with immersion heater,
Power supply: 220-240 volts AC, 50Hz. ;

Power: 480W;

Seamless reservoir with no welds to leak or rust, see- through cover is hinged and removable, and steeply
gabled to accept taller sample

Boiling Chamber, Lid, Cooling Condenser and any other part coming into contact with water / steam made of
quality stainless steel.

Capacity: 25 L

Laboratory Water Bath, Laboratory Water Bath Equipments, Laboratory Water Bath Tools, Laboratory Water
Bath Tool Kits, Laboratory Water Bath Manufacturers, Laboratory Water Bath Suppliers from India, China, Kenya



Laboratory Instrument India