



Email : sales@laboratoryinstrumentindia.com

Product Name :
Six Speed Direct Shift Transmission

Product Code :
LBNY-0005-12000015



Description :

The transmission can be turned easily, the gears changed by hand and compressed air applied in part. The shiftlock is functional. Cutaway are: wet clutch, timing case, transmission case, main shaft, oil cooler, oil filter, oil pump. The transmission is mounted on a table stand. The function of the clutch can be demonstrated by pressing the disks together.

Technical Specification :

Six Speed Direct Shift Transmission



Laboratory Instrument India