



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Six Speed Direct Shift Transmission

**Product Code :**  
LBNY-0005-12000015



**Description :**

The transmission can be turned easily, the gears changed by hand and compressed air applied in part. The shiftlock is functional. Cutaway are: wet clutch, timing case, transmission case, main shaft, oil cooler, oil filter, oil pump. The transmission is mounted on a table stand. The function of the clutch can be demonstrated by pressing the disks together.

**Technical Specification :**

Six Speed Direct Shift Transmission



**Laboratory Instrument India**