



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spring Leg Air Suspension

Product Code :
LBNY-0005-12000016



Description :

The following can be clearly seen on the model: Air suspension withinlet valve for level control and lowering of the chassis, damper with bottomvalve and gas area, solenoid valves for the regulation of the damping strength. In the Daimler Benz S class, four spring legs are used on the front and rear axle asfunction elements of the Airmatic. They are connected with one another via adata bus.

Technical Specification :

Spring Leg Air Suspension



Laboratory Instrument India