




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Direct Shift Transmission

**Product Code :**  
LBNY-0005-12000017



**Description :**

A direct-shift gearbox commonly abbreviated to DSG is an electronically control led dual-clutch multiple-shaft manual gearbox in a transaxle design, without a conventional clutch pedal and with fully automatic or semi-manual control.

Smart Direct Shift Transmission By Getrag:

The following can be seen:

- Synchronization
- Differential
- Duplex clutch
- Parking lock
- Actuators
- Gear wheels
- Mechatronics.

**Technical Specification :**

Direct Shift Transmission

---

# Laboratory Instrument India

