



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Rack And Pinion Steering Gear

**Product Code :**  
LBNY-0005-1200022



**Description :**

Function of thrust block with cup springs:  
The rack is constantly pressed against the pinion, therefore no backlash.

Possible Demonstrations:  
Steering gear ratio  
Movement cycle of steering gear.

**Technical Specification :**

Rack And Pinion Steering Gear



**Laboratory Instrument India**