



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Seven Speed Direct Shift Transmission

**Product Code :**  
LBNY-0005-1200027



**Description :**

The transmission can be turned easily, the gears changed by hand and compressed air applied in part. The shift lock is functional. The transmission is mounted on a table stand.

Cutaway are:

Dry clutch  
Transmission case,  
Oil pump,  
Timing case  
Pressure accumulator  
Main shaft,  
Clutch actuation.

**Technical Specification :**

Seven Speed Direct Shift Transmission



**Laboratory Instrument India**

