



Email : sales@laboratoryinstrumentindia.com

Product Name :
Seven Speed Direct Shift Transmission

Product Code :
LBNY-0005-1200027



Description :

The transmission can be turned easily, the gears changed by hand and compressed air applied in part. The shift lock is functional. The transmission is mounted on a table stand.

Cutaway are:

Dry clutch
Transmission case,
Oil pump,
Timing case
Pressure accumulator
Main shaft,
Clutch actuation.

Technical Specification :

Seven Speed Direct Shift Transmission



Laboratory Instrument India

