



Email : sales@laboratoryinstrumentindia.com

Product Name :
Orbital Shaker

Product Code :
LBNY-0001-510014



Description :

Orbital Shaker

Technical Specification :

Wide range of attachments

Attachment detection

Timer with countdown function

Continuous or touch operation (with standard attachment)

Stable in all speed ranges

Compact, universal small shaker suitable for shaking tasks with all small vessels and micro titer plates, adjustable speed of up to 3,000 rpm.

Type of movement: Orbital movement timer with countdown function.

Shaker speed range: 0-3000rpm

Operating: Button and touch function

Voltage/ frequency input: Voltage: $220 \pm 10\%$ V, 50/60 Hz

Continuous or touch operation: with standard attachment keypad & rotary knob; Timer set up range extended up to 999 min; Low noise.

Orbital Shaker, Orbital Shaker Equipments, Orbital Shaker Tools, Orbital Shaker Tool Kits, Orbital Shaker Manufacturers, Orbital Shaker Suppliers from India, China, Kenya



Laboratory Instrument India