



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Simple Planetary Gear Train

**Product Code :**  
LBNY-0005-12000035



**Description :**

The transmission can be cranked easily. A transmission with superlatives a latest generation transmission. In addition, all the clutches with hydraulic pistons, the crescent pump, the planetary gear trains, the converter with shorting and the hydraulic control device with the solenoid valves. The housing has been cut away to a large extent.

**Technical Specification :**

Simple Planetary Gear Train



**Laboratory Instrument India**