




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Clutch Of The Fourth Generation

**Product Code :**  
LBNY-0005-12000038



**Description :**

The axial piston pumps can also be turned. Everything can be moved the function of the clutch can be well demonstrated by pressing the friction disks together.

The following are easy to recognise:

- Hydraulic pistons
- Controller
- Axial piston pump
- Differential
- Friction disks
- Oil filter and pressure accumulator.

**Technical Specification :**

Clutch Of The Fourth Generation



**Laboratory Instrument India**

