



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cytospin Centrifuge

**Product Code :**  
LBNY-0001-510015



**Description :**

Cytospin Centrifuge, Cytospin Centrifuge Equipments, Cytospin Centrifuge Tools, Cytospin Centrifuge Tool Kits, Cytospin Centrifuge Manufacturers, Cytospin Centrifuge Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara

**Technical Specification :**

Benchtop centrifuge provides economical thin-layer preparations from any liquid matrix, especially hypocellular fluids such as spinal fluid and urine.

- Processes 12 specimens at one time
- Accepts all protocols from Cytospin 1, 2 and 3
- Allows for one-handed opening and closing with a redesigned lid-release mechanism
- Enables viewing of the sealed head through the polycarbonate window during operation
- Allows loading of the sealed head inside a microbiological cabinet for ultimate safety
- Features an auto-locking, plastic outer lid
- Operates only when lid is locked and keeps lid locked during rotation
- Includes disposable sample chambers with caps and a wipe-clean control panel

Technical Specifications:

---

Max. capacity 12 specimens  
Speed range 200 - 2000 min-1  
Rpm or 200-2000 rpm  
Programs 23  
Timer 1 - 99 min  
Anti-unbalance system Audible alarm



**Laboratory Instrument India**