



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electric Motor With Transmission

**Product Code :**  
LBNY-0005-12000041



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The following can be seen:  
Parking lock  
All parts can be turned  
The parking lock is functioning  
Cable harness  
Electric motor  
Differential  
Axle shafts.

**Technical Specification :**

Electric Motor With Transmission



**Laboratory Instrument India**