




Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronically Controlled Shift Valve Body

Product Code :
LBNY-0005-12000045



Description :

An electro-solenoid valve has also been cut away so that its interior can be studied. The Electronically Controlled Shift Valve Body Of An Automatic Transmission gives an insight into the complicated interior of a modern shift- valve body. In the cut-away housing, one sees the control plunger and control springs.

Technical Specification :

Electronically Controlled Shift Valve Body



Laboratory Instrument India