



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Four Speed Automatic Transmission

**Product Code :**  
LBNY-0005-12000048



**Description :**

Electronically controlled.  
This transmission designed for high performance has the following parts cut away:  
All hydraulic couplings  
Torque converter with lockup clutch  
Plan.  
Electronically controlled. The cut-away parts are:  
Planetary gear trains  
Differential and hydraulic control  
Torque converter with lockup clutch  
All hydraulic couplings.  
The transmission can easily be turned.

**Technical Specification :**

Four Speed Automatic Transmission



**Laboratory Instrument India**

