




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Seven Speed Automatic Transmission

**Product Code :**  
LBNY-0005-1200049



**Description :**

The transmission can be turned easily both from the drive and also from the output side. With the 5th generation of automatic transmissions, the following is achieved: high shifting comfort, light weight construction, fuel savings and increased driving enjoyment.

The following are visible:

- All 7 hydraulic couplings
- Range selector can be shifted
- Solenoid valves and many small construction parts
- Parking lock
- Hydraulic control,
- Converter with lock-up clutch
- Ravigneaux planetary gear train and 2 simple planetary geartrains
- Disk springs.

**Technical Specification :**

Seven Speed Automatic Transmission

---

# Laboratory Instrument India

