




Email : sales@laboratoryinstrumentindia.com

Product Name :
Five Speed Automatic Transmission

Product Code :
LBNY-0005-12000055



Description :

Five-speed automatic transmission, electronically controlled
Cut-outs are:
Torque converter with lock up clutch
Shift valve housing with all valves and pistons
Front and rear housing
Crescent-type pump and brake band
All couplings
Internally-gear wheels for planetary gear sets
Actuating piston for brake band.

Technical Specification :

Five Speed Automatic Transmission



Laboratory Instrument India