



Email : sales@laboratoryinstrumentindia.com

Product Name :
ELISA Plate Incubator

Product Code :
LBNY-0001-510017



Description :

ELISA Plate Incubator

Technical Specification :

A distinctive feature of PST plate thermo shaker series is the patented Two-Side Microplates Heating, which allows to achieve full correspondence of the set and actual temperature in the microplate wells.

Position to accommodate at least 2 plates

Micro plate shaker / incubator suitable for all standard depth 96-well plates

Heated lid and base cover the plate entirely

Continuous or timed operation, with alarm buzzer and automatic switch-off

Functional temperature range from -40°C to 120°C

Temperature range: ambient plus 5 to 60?

Temperature stability: approx. 0.1?, uniformity approx. 0.2 ?

Shaking speed: 250 to 1800 rpm, adjustable in steps of 10 rpm

Orbit: 2.mm

LCD displays time set and elapsed, temperature set and actual

Shaking speed: 250 to 1800 rpm, adjustable in steps of 10 rpm

Orbit, 2.mm

Power requirements: 220 V / 50 Hz, with voltage surge protection

ELISA Plate Incubator, ELISA Plate Incubator Equipments, ELISA Plate Incubator Tools, ELISA Plate Incubator Tool Kits, ELISA Plate Incubator Manufacturers, ELISA Plate Incubator Suppliers from India, China, Kenya



Laboratory Instrument India