




Email : sales@laboratoryinstrumentindia.com

Product Name :
Side Valve Four Stroke Engine With Timing Wheel

Product Code :
LBNY-0005-12000066



Description :

The dial gauge can be used to set the upper dead centre exactly on the timing wheel, and thus determine all other piston positions. The engine can be turned easily.

Cutaways are:

Crankcase

Air filter

Valve timing

Cylinder

Cylinder head

Tank

Carburetor and ignition.

Technical Specification :

Side Valve Four Stroke Engine With Timing Wheel



Laboratory Instrument India

