




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Side Valve Four Stroke Engine With Timing Wheel

**Product Code :**  
LBNY-0005-12000066



**Description :**

The dial gauge can be used to set the upper dead centre exactly on the timing wheel, and thus determine all other piston positions. The engine can be turned easily.

Cutaways are:

Crankcase

Air filter

Valve timing

Cylinder

Cylinder head

Tank

Carburetor and ignition.

**Technical Specification :**

Side Valve Four Stroke Engine With Timing Wheel



**Laboratory Instrument India**

