



Email : sales@laboratoryinstrumentindia.com

Product Name :
Screw Compressor With Charge Air Cooling

Product Code :
LBNY-0005-12000067



www.laboratoryinstrumentindia.com

Description :

MG Screw compressor with charge-air cooling

Easily visible:

Two screw spindles, driven via a magnetic clutch and gear wheels.

At the bottom of the charger, is the cutaway charge-aircooler with water channels.

The charger runs at 12,000 RPM.

Technical Specification :

Screw Compressor With Charge Air Cooling



Laboratory Instrument India