



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Screw Compressor With Charge Air Cooling

**Product Code :**  
LBNY-0005-12000067



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

MG Screw compressor with charge-air cooling

Easily visible:

Two screw spindles, driven via a magnetic clutch and gear wheels.

At the bottom of the charger, is the cutaway charge-aircooler with water channels.

The charger runs at 12,000 RPM.

**Technical Specification :**

Screw Compressor With Charge Air Cooling



**Laboratory Instrument India**