



Email : sales@laboratoryinstrumentindia.com

Product Name :
Micro Plate ELISA Washer

Product Code :
LBNY-0001-510018



Description :

Micro Plate ELISA Washer

Technical Specification :

A washer for ELISA microplates designed to ensure thorough washing of reagents between ELISA steps. This equipment is used for many types of ELISA for clinical laboratory applications.

Washer with 8 channel, open system, automatic washer,

Automatic rinse & prime Programme

User definable protocols up to minimum 20

Wash parameters include: 16-character assay name, number of cycles, wash volume, flow rate and variable soak times

Dispense only and aspirate only modes for reagent addition and removal

Built-in multi-speed shaker for improved CVs and reduced assay backgrounds

Crosswise aspiration/double aspiration of flat bottom micro-plates for reduced residual liquid

Bottom wash mode for rapid dilution of reagent

Built-in vacuum & pressure pump assembly

Bottles for waste rinse and wash

Accommodate flat, U or V-shaped bottom plates

Wash cycles between 1-10 wash cycles

Dispensing volumes from 25- to 3000ul

Soak time: 1-600 seconds

Residual well volume should be Fluid flow rate in 150 to 1000ul well/sec to accommodate cellular assays

Spill-over protection and electronics isolated from fluids

Automatic buffer switching

Flip out aerosol cover or similar

All standard accessories, consumables and parts required to operate the equipment, including all standard tools and cleaning and lubrication materials including items not specified above

Operating Temperature: +10 °C to + 30°C

Relative humidity : Electrical Power Supply: 220VAC +10%

Micro Plate ELISA Washer, Micro Plate ELISA Washer Equipments, Micro Plate ELISA Washer Tools, Micro Plate ELISA Washer Tool Kits, Micro Plate ELISA Washer Manufacturers, Micro Plate ELISA Washer Suppliers from India, China, Kenya



Laboratory Instrument India