




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pump With Slide Valve Control

**Product Code :**  
LBNY-0005-12000072



**Description :**

The pump can be connected to a 24 Volt transformer.

Key 1: This moves the control bar. Volume control by operating magnet.

Key 2: This rotates a shaft which moves the slide valve. Injection control by operating magnet.

**Further Functions:**

Control travel sensor,

Pressure valves,

Feed pump,

Emergency-stop valve,

Position sensor and pump pistons.

**Technical Specification :**

Pump With Slide Valve Control



**Laboratory Instrument India**

