




Email : sales@laboratoryinstrumentindia.com

Product Name :
Pump With Slide Valve Control

Product Code :
LBNY-0005-12000072



Description :

The pump can be connected to a 24 Volt transformer.

Key 1: This moves the control bar. Volume control by operating magnet.

Key 2: This rotates a shaft which moves the slide valve. Injection control by operating magnet.

Further Functions:

Control travel sensor,

Pressure valves,

Feed pump,

Emergency-stop valve,

Position sensor and pump pistons.

Technical Specification :

Pump With Slide Valve Control



Laboratory Instrument India

