




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
A Class Common Rail Engine

**Product Code :**  
LBNY-0005-12000074



**Description :**

The engine is driven by a 220V geared motor, all the assemblies also run. All actors and sensors are available and connected to the cable tree. The start can be engaged, the transmission shifted. The following components are cut: cylinder head (with 16 valves), oil injection cooling for pistons, hydraulic valve tappet, oil pump, chain case, generator (internally ventilated), cylinder block, crankcase, oil pan, one piston with cooling duct, starter, exhaust gas recirculation valve, turbocharger, high-pressure pump, one injector, intake manifold, air filter box, feed pump, rail manifold, thermostat, valve cover, air-mass sensor, oil filter with heat exchanger, modern single-vane vacuum pump, 5-gear transmission and differential.

**Technical Specification :**

A Class Common Rail Engine



**Laboratory Instrument India**